



SoClean CPAP cleaner

How does it sanitize your CPAP equipment: Activated oxygen. This is the same technology that is used on produce and drinking water. It disinfects your equipment and naturally breaks down to normal oxygen within two hours. Any excess activated oxygen goes through a filter which converts it back to regular oxygen before release.

Recommended cleaning time: The manufacturer recommends a cleaning cycle of around 7 minutes. After the cleaning cycle is finished, the equipment must remain in the machine for two hours before being worn.

What should I not place in the cleaner: Nylons and natural rubbers (which are not typically used in reservoirs, hoses, masks or head gear). This unit is designed to accept the hose, your mask or nasal pillows, and headgear.

Will my PAP equipment be wet from the SoClean: No -the SoClean disinfects with activated oxygen. No fluids or water are used in this process.

Maintenance: It is recommended that you change the interior filter and check valve system every 6 months.

Where to purchase: www.soclean.com, various CPAP supply websites or by visiting a participating DME. (Search function can be found on SoClean website)

Special Notes: Depending on the brand of CPAP that you are using, you may need to purchase an adapter for your SoClean system. The adapter is included at no cost if you are purchasing the unit from their website. To find out if your machine requires an adapter, [click here](#)

Price: &299.00 – not covered by medical insurance

3B Lumin CPAP Cleaner

How does it sanitize your CPAP equipment: Uses a specialized and high powered UVC light. This is the same type of light source used in hospital operating rooms and sterile clean rooms.

Recommended cleaning time: 5 minutes

What should I not place in the cleaner: According to their user manual, CPAP hoses should not be placed in the cleaner because the optics of the UVC make reaching the interior of the hose difficult. The unit is designed to disinfect CPAP masks, headgear and water chambers.

Will my PAP equipment be wet from the Lumin CPAP cleaner: No -the Lumin uses UV light, not fluids or water.

Maintenance: The Lumin bulb is designed to exceed the life of the device so it should not need changed. If it does become necessary, the bulb can be ordered from the manufacturer and easily replaced by you.

Where to purchase: 3blumin.com or various CPAP supply websites.

Special Notes: Can also be used to sanitize other household items.

Price: \$249.00 – not covered by medical insurance